

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L1V00103	Lot #:	L180301
Description:	Rabbit anti-25-OH-Vitamin D3-BSA (Calcifediol-BSA) Polyclonal Antibody		
Specificity:	Recognizes an epitope on the 25-OH-Vitamin D3-BSA (Calcifediol-BSA)		
Immunogen:	25-OH-Vitamin D3-BSA (Calcifediol-BSA)		
Host Animal:	Rabbit		
Source:	Serum		
Format:	Purified, Liquid		
Concentration:	3.3mg/ml		
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Applications:	ELISA		
Storage:	Store at -20°C for long term storage. Avoid repeated freeze-thaw cycles.		