

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

ProductName:	C-Peptide Antibody	Catalog#:	L3I00603
Description:	PcAb to C-Peptide antibody Rabbit Anti Human C-Peptide Antibody,		
Specificity:	Recognizes an epitope on the C-Peptide Antigen.		
Clone:	None		
Host Animal:	Rabbit		
Isotype:	IgG		
Source:	serum		
Immunogen:	Human C-Peptide Antigen		
Format:	Purified, Liquid		
Purification:	≥90% pure (SDS-PAGE). Protein A chromatography		
Concentration & Lot.:	2.45mg/ml (OD280nm) (Lot. specific)		
Buffer:	PBS, pH 7.4 ± 0.4		
Preservative:	None		
Applications:	Suitable for use in ELISA and Western blot. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelflife:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		