



SPECIFICATION SHEET

Product Name:	Rabbit anti Flu B	Catalog #:	L8H00201
Description:	Rabbit Antibody to Recombinant Human Flu B Antigen		
Host Animal:	Rabbit		
Specificity:	Recombinant Human Flu B Antigen		
Format:	Purified, Liquid		
Concentration&Lot:	11mg/ml	Lot :	L190401
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Applications:	Suitable for use in ELISA, IFA and Western blot. Also suitable for conjugation purposes. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.