

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

ProductName:	Rabbit anti Folic acid (Folate)-BSA	Catalog #:	L6B01103
Synonyms:	Rabbit anti-BSA conjugated FA		
Specificity:	Recognizes an epitope on the Folic acid (Folate)-BSA		
Immunogen:	Folic acid (Folate)-BSA		
Host Animal:	Rabbit		
Source:	Serum		
Format:	Affinity Purified, Liquid		
Concentration:	5.4mg/ml	Lot:	L181101
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		

Description: Folate, distinct forms of which are known as folic acid, folacin, and vitamin B₉, is one of the B vitamins. Synthetic vitamin B₉. Tetrahydrofolate precursor. Promotes DNA synthesis and cell division. Shows anti-inflammatory and chemopreventive effects *in vivo*.

Storage&Shelflife: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.