



SPECIFICATION SHEET

Product Name:	Rabbit anti Hp PcAb	Catalog #:	L1H01202
Description:	Rabbit Antibody to Helicobacter pylori Polyclonal Antibody		
Specificity:	Helicobacter pylori		
Host Animal:	Rabbit		
Immunogen:	Helicobacter pylori		
Format:	Purified, Liquid		
Purification:	>90% pure. Sodium sulfate precipitation and ion-exchange chromatography		
Concentration&Lot:	12mg/ml,	LOT:	L181101
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Applications:	Suitable for use in ELISA, IFA and Western blot. Also suitable for conjugation purposes. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.