



SPECIFICATION SHEET

Product Name:	Recombinant human keratin K19	Catalog #:	L2C00802
Description:	Recombinant human keratin K19 (formerly also designated cytokeratin 19)		
Molecular Weight:	44 098 (calculated from sequence) 40 000 (determined by SDS page)		
Source:	E. coli		
Format:	Purified, Liquid		
Purification:	> 90% (determined by SDS gelelectrophoresis)		
Concentration:	Specific protein concentration: 2mg/ml (CanAg EIA)		
Lot :	L171201		
Buffer:	Tris-HCL, pH 7.4 ± 0.2		
Preservative:	0.2%NaN ₃		
Application:	Immunoassays Immunization.		
Storage:	Store at -20°C. Avoid multiple freeze/thaw cycles.		

This product contains sodium azide, which has been classified as Xn (Harmful) in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.