
SPECIFICATION SHEET

Product Name:	Re combinant cTnl protein	Catalog #:	L4C00102
Description:	Human Troponin I-Cardiac >90% Recombinant, about 30KD		
Source:	Expressed in E. coli		
Format:	Purified		
Purification:	Purified, >90% pure by SDS-PAGE		
Concentration&Lot:	1.1 mg/ml (Roche ECLIA)	Lot :	L20200901
Buffer:	10mmPBS, PH 7.4, 0.01%SKL		
Preservative:	None		
Applications:	Suitable for use as a standard in immunoassay, immunogen for antisera production, and tracer for iodination. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.