

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	HSV-II Ag (Gg)	Catalog #:	L1Y00501
Description:	Recombinant HSV-II (Gg) Antigen Grade		
Source:	Recombinant, E. coli		
Format:	Purified, Protein solution in PBS		
Purification:	Purity >90%		
Concentration& Lot.:	1mg/ml	Lot:	L190401
Buffer:	10Mm PBS, 0.5Mm EDTA pH 7.4± 0.02,		
Preservative:	None		
Applications:	Immune colloidal gold immunochromatography		
Storage& Shelf life:	Upon receipt, store at -20°C. Avoid multiple freeze/thaw cycles.		
Warnings:	For research only		