

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Description: Recombinant human MIS/AMH protein,N-terminal .

Catalog #: L2F00603

Description: Recombinant human anti-Mullerian hormone, AMH

Source: E. Coli **Species:** Homo sapiens

Molecular weight: 77.7kd (with Tag protein)

Concentration: 1.54mg/ml **Lot #:** L181201

Purity: >90% by SDS-PAGE

Formulation: Liquid

Application: ELISA; Other applications have not been tested.
The optimal dilutions should be determined by end user.

Storage buffer: 20mMTrisHCl,PH=8.0, 1mm EDTA,5%Trehalose

Storage & Expiration: Upon receipt, aliquot and store at <-20°C
for long term.Avoid repeated freeze and thaw cycles.

海南省迈迪卡门德商贸有限责任公司
地址: 海南省海口市秀英区秀英大道 18 号宝源花园 6-A-502
邮编: 572000

Hainan Medikamende Trading Co., Ltd.
Adress: 6-A-502, Baoyuanhuayuan,
Xiuyingdadao-No.18, Xiuyingqu,
Haikou, 572000, Hainan, P.R. China

电话: +86-898-31320087 传真: +86-898-31320087 Phone: +86-898-31320087 Fax: +86-898-31320087
网址: www.medical-ivd.com 邮箱: info@medical-ivd.com Website: www.medical-ivd.com Emai: linfo@medical-ivd.com

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

Background:

Müllerian inhibiting substance (MIS), also named anti-Müllerian hormone (AMH), is a tissue-specific TGF-beta superfamily growth factor. Its expression is restricted to the Sertoli cells of fetal and postnatal testis, and to the granulosa cells of postnatal ovary.

海南省迈迪卡门德商贸有限责任公司
地址: 海南省海口市秀英区秀英大道 18 号宝源花园 6-A-502
邮编: 572000

Hainan Medikamende Trading Co., Ltd.
Adress: 6-A-502, Baoyuanhuayuan,
Xiuyingdadao-No.18, Xiuyingqu,
Haikou, 572000, Hainan, P.R. China

电话: +86-898-31320087 传真: +86-898-31320087 Phone: +86-898-31320087 Fax: +86-898-31320087
网址: www.medical-ivd.com 邮箱: info@medical-ivd.com Website: www.medical-ivd.com Emai: linfo@medical-ivd.com