

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikamende Trading Co., Ltd.

SPECIFICATION SHEET

Description: Recombinant human anti-Mullerian hormone Teminal N and C Complex

Catalog #: L2F00604

Description: Coupling Recombinant human anti-Mullerian hormone Teminal N and C Complex

Source: E. Coli **Species:** Homo sapiens

Concentration: 0.2mg/ml (Activation) **Lot #:** L190101

Purity: >90% by SDS-PAGE

Formulation: Liquid

Application: ELISA; Other applications have not been tested.
The optimal dilutions should be determined by end user.

Storage buffer: 10mM CBS PH=9.0, 1mm EDTA,5%Trehalose

Storage & Expiration: Upon receipt, aliquot and store at <-20°C
for long term.Avoid repeated freeze and thaw cycles.

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

Background:

Müllerian inhibiting substance (MIS), also named anti-Müllerian hormone (AMH), is a tissue-specific TGF-beta superfamily growth factor. Its expression is restricted to the Sertoli cells of fetal and postnatal testis, and to the granulosa cells of postnatal ovary.