

## SPECIFICATION SHEET

<b>Product Name:</b>	B19 VP1+2	<b>Catalog #:</b>	L4H002
<b>Description:</b>	Recombinant human B19 VP1+VP2 Antigen (B19 VP1+VP2)		
<b>Molecular Weight:</b>	33KD (calculated from sequence) 33KD (determined by SDS gelelectrophoresis)		
<b>Source:</b>	E. coli		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	> 90% (determined by SDS page)		
<b>Concentration&amp;Lot:</b>	2.5 mg/ml	<b>Lot :</b>	L170701
<b>Buffer:</b>	0.01M PBS pH 7.4 ± 0.2 0.5mM EDTA		
<b>Preservative:</b>	None		
<b>Application:</b>	Immunoassays as EIA, IFA or Immunization.		
<b>Storage&amp;Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		