

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

ProductName: GAD **Catalog #:** L4I00302

Description: Recombinant human Glutamic Acid Decarboxylase(GAD65) Antigen

Molecular Weight: 80kd (calculated from sequence)
80KD(determined by SDS gelelectrophoresis)

Source: E. coli

Format: Purified, Liquid

Purification: > 85% (determined by SDS gelelectrophoresis)

Concentration& Lot.: 1.5mg/ml (OD280) **Lot:** L180801

Buffer: Tris-HCL, pH 7.4 ± 0.2

Preservative: 0.02%NaN3

Application: Immunoassays Immunization.

Storage&Shelflife:

Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.