

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Product Name:** GAD65 **Catalog #:** L4I00304

**Description:** Recombinant human Glutamic Acid Decarboxylase(GAD2/GAD65)

**Molecular Weight:** 65kd (calculated from sequence)  
70KD(determined by SDS gelelectrophoresis)

**Source:** Pichia pastoris

**Format:** Purified, Liquid

**Purification:** > 85% (determined by SDS gelelectrophoresis)

**Concentration& Lot.:** 1mg/ml **Lot:** L151201

**Buffer:** 50mM Tris, 5% glycerol pH:8.0

**Preservative:** None

**Application:** Immunoassays Immunization.

**Storage&Shelf life:** Store at -20°C for 3 years. Avoid multiple