

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

ProductName: IA2 **Catalog#:** L4I00501

Description: Recombinant human insulinoma associated protein 2 (IA2)

MW: 60KD

Source: E. coli

Format: Purified, Liquid

Purification: > 85% (determined by SDS gelelectrophoresis)

Concentration& Lot.: 1.5 mg/ml **Lot.:** L180401

Buffer: 10mmPBS,1mmEDTA,pH7.4

Preservative: None

Storage&Shelflife: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.