

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Description:** Recombinant human MIS/AMH protein. **Catalog #:** L2F00602

**Description:** Recombinant human anti-Mullerian hormone, AMH

**Source:** E. Coli

**Molecular weight:** 15 KD

**Concentration:** 2 mg/ml **Lot #:** L170401

**Purity:** >90% by SDS-PAGE

**Formulation:** Liquid

**Application:** ELISA; Other applications have not been tested.  
The optimal dilutions should be determined by end user.

**Storage buffer:** 50mMTris、PH=8.0, 1mm EDTA,1mMDTT,0.05% SKL

**Storage & Expiration:** Upon receipt, aliquot and store at <-20°C

for long term.Avoid repeated freeze and thaw cycles.

#### Background:

Müllerian inhibiting substance (MIS), also named anti-Müllerian hormone (AMH), is a tissue-specific TGF-beta superfamily growth factor. Its expression is restricted to the Sertoli cells of fetal and postnatal testis, and to the granulosa cells of postnatal ovary.