



## SPECIFICATION SHEET

**Product Name:** Recombinant Human NuMA Protein(nmp22)

**Catalog #:** L2C03701

**SYNONYMS:** centrophilin stabilizes mitotic spindle in mitotic cells、NMP-22、nuclear mitotic apparatus protein 1、NuMA protein、NuMA、NUMA1、SP-H Antigen、structural nuclear protein

**Source:** E. Coli **Species:** Homo sapiens

**Molecular weight:** 50kd (with Tag protein)

**Concentration:** 2.18mg/ml

**Lot #:** L20211201

**Purity :** >90% by SDS-PAGE

**Formulation :** The Protein solution contains 10mMPBS,PH=7.4, 1mm EDTA, 1mM DTT, 5% Trehalose.

**Application:** ELISA; Positive Control; Immunogen; SDS-PAGE; WB. Other applications have not been tested .

**Storage & Expiration:** Upon receipt, aliquot and store at <-20℃ .for long term.Avoid repeated freeze and thaw cycles.

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**