

SPECIFICATION SHEET

Description:	Recombinant human SCCA Antigen Catalog #: L2C01501
Molecular Weight:	48kd (calculated from sequence) 48 KD(determined by SDS gelelectrophoresis)
Source:	E. coli
Format:	Purified, Liquid
Concentration:	Specific protein concentration: 1.8mg/ml (CanAg EIA)
Lot:	L190302
Buffer:	20mm Tris-HCL, pH 7.0 + 0.2, 0.5mm EDTA
Preservative:	0.02%NaN3
Storage&Shelf life:	Store at -20°C to 2021-7. Avoid multiple freeze/thaw cycles

This product contains sodium azide, which has been classified as Xn (Harmful) in European Directive 67/548/EEC in the concentration range of 0.1 – 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.