



SPECIFICATION SHEET

Description: Recombinant human SCCA2 Antigen **Catalog #:**L2C01503

Molecular Weight: 48kd (calculated from sequence)
48 KD(determined by SDS gelelectrophoresis)

Source: E. coli

Format: Purified, Liquid

Purification: > 90% (determined by SDS gelelectrophoresis)

Concentration: Specific protein concentration: 0.8 mg/ml (CanAg EIA)

Lot: L150701

Buffer: Tris-HCL, pH 7.4 ± 0.2

Preservative: 0.02%NaN₃

Application: Immunoassays Immunization.

Storage: Store at -20°C. Avoid multiple freeze/thaw cycles.