



SPECIFICATION SHEET

Catalog #: L4I01101

Description: Recombinant RBP-4 (Retinol-Binding Protein 4,RBP 4)

SYNONYMS: Retinol Binding Protein 4, RBP-4, RBP4, Plasma retinol-binding protein, PRBP, RBP.

Source: E. Coli **Species:** Homo sapiens

Molecular weight: 21kd (with His Tag protein)

Concentration: 1mg/ml

Lot #: L170801

Purity : >95% by SDS-PAGE

Formulation : The Protein solution contains 20mM saline sodium citrate pH3.5, 1mm EDTA,5%Trehalose

Preservative: None

Application: ELISA; Positive Control; Immunogen; SDS-PAGE; WB. Other applications have not been tested . If bio-activity of the protein is needed, please check active protein.

Storage & Expiration: Upon receipt, aliquot and store at <-20°C .for long term.Avoid repeated freeze and thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.