



SPECIFICATION SHEET

Product Name:	TB CFP-10	Catalog #:	L1H01304
Description:	Recombinant TB CFP-10 antigen, Antigen Grade		
Source:	Recombinant		
Format:	Purified, Protein solution in PB		
Purification:	Purity >90%		
Concentration&Lot:	2.5mg/ml	Lot:	L180501
Buffer:	10Mm PB, pH 7.4 ± 0.02		
Preservative:	None		
Applications:	EIA,CLIA,IP,IHC,WB		
Storage&Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
Warnings:	For research only		