

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**Catalog #:** L2K002

**Lot #: L190101**

**Description:** S-(5'-Adenosyl)-L-homocysteine is S - adenosine methionine dependence methyltransferase reaction in the process of competitive inhibition. It is also known as 5'-Deoxy-S-adenosyl-L-homocysteine AdoHcy ;  
S-(5'-Deoxyadenosine-5')-L-homocysteine .

**Appearance:** White or almost white crystalline powder; odorless; tasteless.

**CAS:** 979-92-0

**Formula:**  $C_{14}H_{20}N_6O_5S$

**MW:** 384.41

**Form:**

Crystalline, Redissolve by 1M HCl:soluble 19.6-20.40mg/ml ,  
clear to slightly hazy,colorless to faintly yellow

**Contaminants:** not detectable

**Storage:** -20°C or lower, for long term storage.