

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name: Tox-Ag P30 **Catalog #:** L1Y00301

Description: Recombinant T. gondii P30 Ag (P30,SAG1)

Source: E.Coli, 30KD
Strain: RH from ATCC

Format: Lysate, Liquid

Purification: >85% (SDS-PAGE)

Concentration&Lot.: 2mg/ml (Lot. L161001)

Buffer: PBS, 0.15M NaCl, pH 7.4 \pm 0.2

Preservative: None

Applications: Specific methodologies have not been tested using this product.

Storage&Shelf life:

Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.