

---

SPECIFICATION SHEET

<b>Product Name:</b>	McAb-HIV-1P24	<b>Catalog #:</b>	L2H00503
<b>Description:</b>	MAb to HIV-1 p24 Monoclonal Antibody to Human Immunodeficiency Virus Type 1 (HIV-1), p24 specific		
<b>Specificity:</b>	HIV-1 p24 (Recombinant antigen)		
<b>Clone:</b>	MLP3E12		
<b>Host Animal:</b>	Mouse		
<b>Isotype:</b>	IgG		
<b>Source:</b>	Ascites		
<b>Immunogen:</b>	Recombinant Antigen		
<b>Format:</b>	Purified, Liquid		
<b>Concentration:</b>	1.8mg/ml (OD280nm)	<b>Lot:</b>	L190401
<b>Buffer:</b>	0.01M PBS, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	Suitable as a detection antibody in ELISA. Each laboratory should determine an optimum working titer for use in its particular application.		
<b>Storage:</b>	Store at -20°C. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.