

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

**Description:** L-T<sub>3</sub> 3,3', 5-Triiodo-L-thyronine Sodium Salt

**Catalog #:** L2G01001

**Molecular Formula:** C<sub>15</sub>H<sub>11</sub>I<sub>3</sub>NNaO<sub>4</sub>

**MW:** 672.96

**Purity:** 95%

**Lot:** L160601

**Forma:** Powder

**Storage:** Store at 2 ~ 8°C, Against light or wet.

**Dissolved:** Dissolve one by one drops of NaOH 0.1mol/L, and make it dissolve completely. The solution is diluted to the appropriate concentration with 50% propylene glycol solution.(Note: this method is only applicable to the calibration of the product, for reference. Other uses of the dissolution method please refer to the relevant information.)