



SPECIFICATION SHEET

Catalog #: L1K00102 **Lot #:** L160301

Description: McAb to Cystatin C antibody
Mouse Anti Human Cystatin C Monoclonal Antibody

Specificity: Recognizes an epitope on Cystatin C Antigen

Clone: CY002

Host Animal: Mouse

Source: Ascites

Immunogen: Recombinant Cystatin C Antigen

Format: Purified, Liquid

Purification: Protein A chromatography

Concentration: 3.5mg/ml

Buffer: 0.01M PBS, 0.15M NaCl, pH 7.4

Preservative: None

Storage&Shelf life: Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.