

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

ProductName:	McAb to S100-A8	Catalog #:	L6B00402
Description:	Mouset anti Recombinant Human S100 Calcium Binding Protein A8		
Specificity:	Recognizes an epitope on Human S100 Calcium Binding Protein A8		
Immunogen:	Recombinant Human S100 Calcium Binding Protein A8		
Host Animal:	Mouse		
Source:	Serum		
Format:	Purified, Liquid		
Concentration:	5.8mg/ml	LOT:	L180801
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Storage&Shelflife:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		