

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Product Name:</b>	McAb to AMH	<b>Catalog #:</b>	L1F00803				
<b>Description:</b>	Monoclonal Antibody to Human anti-Mullerian hormone, AMH						
<b>Clone:</b>	AM003						
<b>Host Animal:</b>	Mouse						
<b>Isotype:</b>	IgG						
<b>Source:</b>	Ascites						
<b>Immunogen:</b>	Recombinant Human anti-Mullerian hormone, AMH						
<b>Format:</b>	Purified, Liquid						
<b>Concentration&amp;Lot:</b>	4.66 mg/ml (OD280nm)	<b>Lot #:</b>	L171201				
<b>Buffer:</b>	10mM PBS pH 7.4 ± 0.2						
<b>Preservativ:</b>	None						
<b>Applications:</b>	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1F00801</td><td>L1F00803</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1F00801	L1F00803
<u>Capture</u>	<u>Detection</u>						
L1F00801	L1F00803						
<b>Storage:</b>	Store at -20°C for at least 3 years.						