

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Catalog #:</b>	L1G00801	<b>Lot #:</b> L180801				
<b>Description:</b>	MAb to FSH Monoclonal Antibody to Human <b>Follicle Stimulating Hormone (FSH)</b> alpha subunit					
<b>Clone:</b>	FS001					
<b>Host Animal:</b>	Mouse					
<b>Isotype:</b>	IgG <sub>1</sub> ,					
<b>Source:</b>	Ascites					
<b>Immunogen:</b>	Purified natural FSH alpha					
<b>Format:</b>	Purified, Liquid					
<b>Concentration:</b>	1.7 mg/ml					
<b>Buffer:</b>	10mM PBS pH 7.4 ± 0.2					
<b>Preservative:</b>	None					
<b>Applications:</b>	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00801</td><td>L1G00802</td></tr></table>		<u>Capture</u>	<u>Detection</u>	L1G00801	L1G00802
<u>Capture</u>	<u>Detection</u>					
L1G00801	L1G00802					
<b>Storage:</b>	Short term (up to 2 Weeks) store at 2–8°C. Long term, store at <-20°C.					