

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

| | | | | | | |
|-----------------------|---|-----------------------|----------------|------------------|----------|----------|
| Catalog #: | L1G00802 | Lot #: L180801 | | | | |
| Description: | MAb to FSH Monoclonal Antibody to Human Follicle Stimulating Hormone (FSH) beta subunit | | | | | |
| Clone: | FS002 | | | | | |
| Host Animal: | Mouse | | | | | |
| Isotype: | IgG | | | | | |
| Source: | Ascites | | | | | |
| Immunogen: | Purified natural FSH beta | | | | | |
| Format: | Purified, Liquid | | | | | |
| Concentration: | 2 mg/ml | | | | | |
| Buffer: | 10mM PBS pH 7.4 ± 0.2 | | | | | |
| Preservative: | None | | | | | |
| Applications: | Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00801</td><td>L1G00802</td></tr></table> | | <u>Capture</u> | <u>Detection</u> | L1G00801 | L1G00802 |
| <u>Capture</u> | <u>Detection</u> | | | | | |
| L1G00801 | L1G00802 | | | | | |
| Storage: | Store at <-20°C for at least 3 years. | | | | | |