



SPECIFICATION SHEET

**Product Name:** McAb to HBcAg                      **Catalog #:** L1H00601

**Description:** Monoclonal Antibody to Hepatitis B c Antigen (HBcAg)

**Specificity:** HBcAg.

**Clone:** LT05

**Host Animal:** Mouse. Hybridization of Sp2/0 myeloma

**Isotype:** IgG cells with spleen cells from BALB/c mice.

**Source:** Ascites

**Immunogen:** HBcAg

**Format:** Purified, Liquid

**Purification:** >90% pure. Protein A chromatography

**Concentration:** 2.5mg/ml      **Lot.:** H141001

**Buffer:** 10MM PBS, pH 7.4

**Preservative:** None

**Applications:** Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.

**Storage& Shelf life:** Store at -20° C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.