



SPECIFICATION SHEET

Product Name:	McAb to HBeAg	Catalog #:	L1H00401				
Description:	Monoclonal Antibody to Hepatitis B e Antigen (HBeAg)						
Specificity:	HBeAg.						
Clone:	N/A						
Host Animal:	Mouse. Hybridization of Sp2/0 myeloma						
Isotype:	IgG cells with spleen cells from BALB/c mice.						
Source:	Ascites						
Immunogen:	Native HBeAg						
Format:	Purified, Liquid						
Purification:	>90% pure. Protein A chromatography						
Concentration:	3.5mg/ml	Lot. :	L180501				
Buffer:	PBS, pH 7.4						
Preservative:	None						
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1H00401</td><td>L1H00402</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1H00401	L1H00402
<u>Capture</u>	<u>Detection</u>						
L1H00401	L1H00402						

Storage& Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende

Trading Co., Ltd.

地址：海南省海口市秀英区秀英大道18号宝源花园6-A-502

Address: 6-A-502, Baoyuanhuayuan,

Xiuyingdadao-No. 18,

邮编：572000