

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L1G00501	Lot #: L190301				
Description:	MAb to alfa-LH,HCG,FSH Monoclonal Antibody to Human alfa submit of Luteinizing Hormone.					
Specificity:	Human alfa-LH Antigen					
Clone:	LH001					
Host Animal:	Mouse					
Source:	Ascites					
Immunogen:	Purified natural LH, alfa chain					
Format:	Purified, Liquid					
Purification:	≥90% pure (SDS-PAGE).					
Concentration:	5mg/ml					
Buffer:	10mM PB, pH 7.4 ± 0.2 containing 150mM NaCl					
Preservative:	None					
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00501</td><td>L1G00502</td></tr></table>		<u>Capture</u>	<u>Detection</u>	L1G00501	L1G00502
<u>Capture</u>	<u>Detection</u>					
L1G00501	L1G00502					
Storage:	Store at <-20°C for three years.					