

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L1G00502	Lot #: L190301				
Description:	MAb to beta-LH Monoclonal Antibody to Human beta-LH .					
Specificity:	Human beta-LH Antigen					
Clone:	LH002					
Host Animal:	Mouse					
Isotype:	IgG ₁ ,					
Source:	Ascites					
Immunogen:	Purified natual LH					
Format:	Purified, Liquid					
Purification:	≥90% pure (SDS-PAGE). Protein A chromatography.					
Concentration:	5.4 mg/ml					
Buffer:	10mM PB, pH 7.4 ± 0.2 containing 150mM NaCl					
Preservative:	None					
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00501</td><td>L1G00502</td></tr></table>		<u>Capture</u>	<u>Detection</u>	L1G00501	L1G00502
<u>Capture</u>	<u>Detection</u>					
L1G00501	L1G00502					
Storage:	Store at <-20°C for three years.					