

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikamende Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L1G00101	Lot #: L160701				
Description:	McAb to PRL Monoclonal Antibody to prolactin, PRL.					
Specificity:	Human prolactin, PRL					
Clone:	PR001					
Host Animal:	Mouse					
Source:	Ascites					
Immunogen:	Purified Recombinant prolactin					
Format:	Purified, Liquid					
Purification:	≥90% pure (SDS-PAGE).					
Concentration:	2.1mg/ml					
Buffer:	10mM PB, pH 7.4 ± 0.2 containing 150mM NaCl					
Preservative:	None					
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00101</td><td>L1G00102</td></tr></table>		<u>Capture</u>	<u>Detection</u>	L1G00101	L1G00102
<u>Capture</u>	<u>Detection</u>					
L1G00101	L1G00102					
Storage:	Store at <-20°C for three years.					