

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Description:	Mouse anti- transferrin, TRF antibody Catalog#: L1K01001						
Background:	Transferrins are iron-binding blood plasma glycoproteins that control the level of free iron (Fe) in biological fluids. Human transferrin is encoded by the TF gene.						
Specificity:	Human						
Clone:	TRF001						
Host Animal:	Mouse.						
Isotype:	IgG						
Immunogen:	Purified Human Serum Albumin						
Format:	Purified, Liquid						
Purification:	>90% pure. Chromatographically purified and supplied in PBS, pH 7.4.						
Concentration&Lot.:	2.4 mg/ml Lot: L190501						
Buffer:	PBS, pH 7.4						
Applications:	Detection of insulin of two-site enzyme immunoassay. Each laboratory should determine an optimum working titer for use in its particular application. The following pair is recommended for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L3I00701</td><td>L1K00401</td></tr><tr><td>L1K00401</td><td>L1K00402, L1K00403</td></tr></table>	<u>Capture</u>	<u>Detection</u>	L3I00701	L1K00401	L1K00401	L1K00402, L1K00403
<u>Capture</u>	<u>Detection</u>						
L3I00701	L1K00401						
L1K00401	L1K00402, L1K00403						
Storage&Shelflife:	Store at -20°C for 3 years. Avoid multiple freeze/thaw						