

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikamende Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L1G00601	Lot #: L170701				
Description:	MAb to TSH Monoclonal Antibody to Human Thyroid Stimulating Hormone (TSH) beta subunit					
Clone:	TS001					
Host Animal:	Mouse					
Isotype:	IgG ₁ ,					
Source:	Ascites					
Immunogen:	Purified natural TSH beta					
Format:	Purified, Liquid					
Purification:	Ion exchange $\geq 90\%$ pure.					
Concentration:	1.9 mg/ml					
Buffer:	10mM PBS pH 7.4 \pm 0.2, 0.01% NaN ₃					
Preservative:	0.01% NaN ₃					
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pairs for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00601</td><td>L1G00602</td></tr></table>		<u>Capture</u>	<u>Detection</u>	L1G00601	L1G00602
<u>Capture</u>	<u>Detection</u>					
L1G00601	L1G00602					
Storage:	Store at -20°C for three years.					