

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L1G00602	Lot #:	L151101				
Description:	MAb to TSH Monoclonal Antibody to Human Thyroid Stimulating Hormone (TSH), beta subunit						
Clone:	TS002						
Host Animal:	Mouse						
Isotype:	IgG ₁						
Source:	Ascites						
Immunogen:	Purified natural TSH beta						
Format:	Purified, Liquid						
Purification:	Ion exchange $\geq 90\%$ pure.						
Concentration:	3.2 mg/ml						
Buffer:	10mM PBS pH 7.4 \pm 0.2, 0.01% NaN ₃						
Preservative:	0.01% NaN ₃						
Applications:	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Recommended pair for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1G00601</td><td>L1G00602</td></tr></table>			<u>Capture</u>	<u>Detection</u>	L1G00601	L1G00602
<u>Capture</u>	<u>Detection</u>						
L1G00601	L1G00602						
Storage:	Store at -20°C for three years.						