

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

ProductName: G-17-Ag-BSA **Catalog #:** L5B01001

Description: BSA Conjugated Human GASTRIN 17 (G-17)

Format: liquid, G-17 purified by HPLC

Purification: Purified

Concentration& Lot.: 10mg/ml **Lot: L190901**

Buffer: 10Mm PBS, pH 7.5

Preservative: None

Description:

Endogenous peptide, mainly found in the gut, that activates gastric acid secretion, delays gastric emptying and increases blood flow in the gastric mucosa.

Storage&Shelflife: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.