

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

<b>Catalog #:</b>	L2K00201	<b>Lot #:</b>	L180501
<b>Description:</b>	BSA-S-adenosyl-L-homocysteine BSA conjugated SAH		
<b>Type:</b>	Synthetic		
<b>Format:</b>	Purified, Liquid		
<b>Concentration:</b>	7.7mg/ml		
<b>Buffer:</b>	0.01M PBS, pH 7.4		
<b>Preservative:</b>	None		
<b>Applications:</b>	ELISA		
<b>Storage:</b>	Store at -20°C for long term storage. Avoid repeated freeze-thaw cycles.		