

**SPECIFICATION SHEET**

Product Name:	Human CA24-2 Calib. Grade	Catalog #:	L2C00602
Description:	Human Cancer Antigen 24-2 (CA24-2) (Gastrointestinal Cancer) Calibration Grade		
Source:	Human colion carcinoma.		
Format:	Purified, Liquid		
Purification:	Gel filtration and column chromatography.		
Concentration:	450KU/ml (Canag ELISA)	Lot #:	L170201
Buffer:	PBS, pH 7.4		
Preservative:	0.05% Sodium Azide		
Applications:	Immunogen, calibrator in immunoassays and in vitro studies. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.		
Storage:	Upon receipt, store at -20°C.		
Inactivation:	Not applicable.		
Warnings:	All human source materials have tested negative for HIV1, HIV2, HCV antibodies and HBsAg. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		