

# 海南省迈迪卡门德商贸有限责任公司

## Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

---

### SPECIFICATION SHEET

**ProductName:** CMV Angiten **Catalog#:** L1Y00201

**Description:** Cytomegalovirus Ag (P52+PP150)  
Cytomegalovirus M Antigen. (P52+PP150)  
WM: 25KD

**Source:** Recombinant ,E.col

**Format:** Purified, Protein solution in Tris-HCl pH8.0,

**Purification:** Purity >85%

**Concentration&Lot.:** 2mg/ml **Lot.:** L180701

**Buffer:** PBS, pH 7.4 + 0.2 , 0.5mM EDTA

**Preservative:** None

**Applications:** Immune colloidal gold immunochromatography

**Storage&Shelf life:**

Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

**Warnings:** For research only

网址: [www.medical-ivd.com](http://www.medical-ivd.com) 邮箱: [info@medical-ivd.com](mailto:info@medical-ivd.com)

Website: [www.medical-ivd.com](http://www.medical-ivd.com) Emai: [linfo@medical-ivd.com](mailto:linfo@medical-ivd.com)