

海南省迈迪卡门德商贸有限责任公司
Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #: L2G00201 **Lot #: L170501**

Description: Follicle Stimulating Hormone (FSH) intact human

Purity: >85% by SDS-PAGE

Source: Human pituitary glands

Concentration: 1mg/ml (1800U/mg)

Contaminants: not detectable

Buffer: 10mM PBS, 150mM NaCl, PH 7.4

Biohazard: Raw materials tested negative against Syphilis, HCV, HIV 1/2, HBcAg and HBsAg. Use full laboratory safety precautions.

Storage: -20°C or lower, for long term storage.