
SPECIFICATION SHEET

| | | | | | | | |
|-----------------------|---|-------------------|----------|----------------|------------------|----------|----------|
| Product Name : | Goat Anti HBsAg PcAb | Catalog #: | L1H00201 | | | | |
| Description: | Goat Antibody to Human Hepatitis B Surface Antigen (HBsAg) | | | | | | |
| Specificity: | Reacts with both ad and ay antigens. No reaction when tested against normal human serum and single arc when tested against HBsAg elevated human serum by IEP. | | | | | | |
| Host Animal: | Goat | | | | | | |
| Immunogen: | Purified Human Hepatitis B Surface Ag | | | | | | |
| Format: | Purified, Liquid | | | | | | |
| Concentration: | 3.8 mg/ml | Lot. : | L180601 | | | | |
| Buffer: | 10Mm PBS, pH 7.4 | | | | | | |
| Preservative: | None | | | | | | |
| Applications: | Suitable for use in ELISA. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Recommended pair for sandwich immunoassay: <table><tr><td><u>Capture</u></td><td><u>Detection</u></td></tr><tr><td>L1H00202</td><td>L1H00201</td></tr></table> | | | <u>Capture</u> | <u>Detection</u> | L1H00202 | L1H00201 |
| <u>Capture</u> | <u>Detection</u> | | | | | | |
| L1H00202 | L1H00201 | | | | | | |

Storage& Shelf life: Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.