



Hainan Medikamende Trading Co., Ltd.

海南省迈迪卡门德商贸有限责任公司
Hainan Medikamende Trading Co., Ltd.

SPECIFICATION SHEET

Product Name:	Goat anti human IgA	Catalog #:	L1102102
Description:	Goat Antibody to human Immunoglobulin A (IgA).		
Specificity:	Reacts with heavy chains on human IgA and with the light chains on all human immunoglobulins. By Immunoelectrophoresis (IEP) shows bands only in the gamma region against Normal human Serum.		
Host Animal:	Goat		
Immunogen:	Purified human IgA		
Format:	Affinity Purified, Liquid		
Purification:	>90% pure (SDS-PAGE).		
Concentration:	1.55mg/ml	Lot.:	L161201
Buffer:	10mM Sodium Phosphate, 150mM Sodium Chloride, pH 7.2		
Preservative:	None		
Applications:	Specific methodologies have not been tested using this product.		
Storage& Shelf life:	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		

FOR RESEARCH USE ONLY

海南省迈迪卡门德商贸有限责任公司
地址：海南省海口市秀英区秀英大道18号宝源花园6-A-502
邮编：572000
电话：+86-898-31320087，传真：+86-898-31320087
网址：www.medical-ivd.com 邮箱：info@medical-ivd.com

Hainan Medikamende Trading Co., Ltd.
Address: 6-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,
Xiuyingqu, Haikou, 572000, Hainan, P.R. China
Phone: +86-898-31320087 Fax: +86-898-31320087
Website: www.medical-ivd.com Emai: linfo@medical-ivd.com