

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Catalog #:	L1G00302	Lot #:	L180801
Description:	PcAb to Thyroxine (T4) Goat anti Thyroxine (T4) Ponoclonal Antibody		
Specificity:	<u>Cross-reactivity:</u> T3<1%		
Clone:	N/A		
Host Animal:	Goat		
Isotype:	IgG		
Source:	Ascites		
Immunogen:	T4-BSA		
Format:	Purified, Liquid		
Purification:	>90% pure. Chromatographic methods		
Concentration:	1.7mg/ml		
Buffer:	0.01M PBS, pH 7.4		
Preservative:	None		
Applications:	Specific methodologies have not been tested using this product.		
Storage:	store at <-20°C.		