



---

## SPECIFICATION SHEET

<b>Product Name:</b>	Hepatitis B Surface Ag (HBs Ag)	<b>Catalog #:</b>	L1H00301
<b>Description:</b>	Hepatitis B Surface Antigen		
<b>Source:</b>	Human infected plasma		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	Ion exchange chromatography		
<b>Concentration:</b>	1.75 mg/ml	<b>Lot. :</b>	L181201
<b>Buffer:</b>	0.01M PBS, pH 7.4		
<b>Preservative:</b>	0.02% Sodium azide		
<b>Applications:</b>	Specific methodologies have not been tested using this product.		
<b>Storage&amp; Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
<b>Inactivation:</b>	Heated for 10 hours at 60°C before ultracentrifugation.		
<b>Warnings:</b>	All human source materials have tested negative for HIV1, HIV2 and HCV antibodies. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		
<b>For Research Use Only.</b>			