



---

## Specification Sheet

<b>Product Name:</b>	HCV MP Antigen	<b>Catalog # :</b>	L1H00801
<b>Description:</b>	Hepatitis C Virus HCV (Core+NS3+NS4+NS5) Antigen, genotypecore+3+4+5.		
<b>Source:</b>	E. coli		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	95% pure (SDS-PAGE)		
<b>Concentration:</b>	1.6mg/ml (OD280nm).(Lot. specific)		
<b>Buffer:</b>	10Mm PBS PH7.4		
<b>Preservative:</b>	none		
<b>Applications:</b>	Specific methodologies have not been tested using this product.		
<b>Storage&amp; Shelf life:</b>	Store at -20°C for 3 years. Avoid multiple freeze/thaw cycles.		
<b>Warnings:</b>	All human source materials have tested negative for HIV1, HIV2 and HCV antibodies. No test guarantees a product to be non-infectious. Therefore, all material derived from human fluids or tissues should be considered as potentially infectious.		

### For Research Use Only.