

## SPECIFICATION SHEET

**Product Name:** HIV-1 P24                      **Catalog #:** L2H002

**Description:** HIV-1 P24 Recombinant antigen, Antigen Grade

**Source:** Recombinant ,E.coli

**Format:** Purified, Protein solution in Tris-HCl pH8.0,

**Purification:** Purity >90%

**Concentration&Lot:** 2mg/ml                      **Lot:** L181201

**Buffer:** 20Mm Tris-HCl, pH 8.0 ± 0.02,

**Preservative:** None

**Application:** Immunoassays or Immunization.

**Storage&Shelf life:** Store at -20 °C for 3 years. Avoid multiple freeze/thaw cycles.

**Warnings:** For research only