



---

## SPECIFICATION SHEET

<b>Product Name:</b>	Recombinant HP UreaseA	<b>Catalog #:</b>	L1H01203
<b>Description:</b>	Rcombinant Helicobacter pylori (Hp) Urease A Antigen, Antigen Grade		
<b>Molecular Weight:</b>	26.71KD (calculated from sequence) 26.7KD (determined by SDS gelelectrophoresis)		
<b>Source:</b>	E. coli		
<b>Format:</b>	Purified, Liquid		
<b>Purification:</b>	> 85% (determined by SDS gelelectrophoresis)		
<b>Concentration&amp;Lot:</b>	1.29mg/ml,	<b>Lot:</b>	L180401
<b>Buffer:</b>	20Mm Tris-HCl, pH 8.0 ± 0.1		
<b>Preservative:</b>	0.2%NaN3		
<b>Application:</b>	Immunoassays or Immunization.		
<b>Storage&amp;Shelf life:</b>	Store at -20 C for 3 years. Avoid multiple freeze/thaw cycles.		

---

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.

地址：海南省海口市秀英区秀英大道18号宝源花园6-A-502

Address: 6

-A-502, Baoyuanhuayuan, Xiuyingdadao-No.18,

邮编：572000

Xiuyingqu, Haikou, 572000, Hainan, P.R. China

电话：+86-898-31320087，传真：+86-898-31320087

Phone: +86-898-31320087

Fax: +86-898-31320087

网址：www.medical-ivd.com 邮箱：info@medical-ivd.com

Website: www.medical-ivd.com

Emai: linfo@medical-ivd.com