

海南省迈迪卡门德商贸有限责任公司

Hainan Medikamende Trading Co., Ltd.



Hainan Medikament Trading Co., Ltd.

SPECIFICATION SHEET

Product Name: HSV-1 gG protein **Catalog #:** L1Y00401

Description: Herpes Simplex Virus-1, gG protein, 46KD
Herpes Simplex Virus Type 1 (HSV-1), gG protein

Source: Recombinant, E.coli

Format: Purified, Liquid

Purification: Affinity purified, SDS-Page >90%

Concentration& Lot.: 3mg/ml (OD280)(lot. L160301)

Buffer: 10Mm PBS, 2Mm EDTA, pH 7.4± 0.02,

Preservative: None

Applications: Immune colloidal gold immunochromatography

Storage& Shelf life: Upon receipt, store at -20°C.
Avoid multiple freeze/thaw cycles.